



Elabscience®

www.elabscience.com

ELISA Kits For Cancer Research



Cancer Research

There are more than 100 types of cancer, which are usually named for the organs or tissues where the cancers form. In all types of cancer, some of the body's cells begin to divide without stopping and spread into surrounding tissues.

Cancer can start almost anywhere in the body, which is made up of trillions of cells. Cancer cells are able to ignore signals that normally tell cells to stop dividing or that begin a process known as programmed cell death, or apoptosis, which the body uses to get rid of unneeded cells. Cancer cells are also often able to evade the immune system, a network of organs, tissues, and specialized cells that protects the body from infections and other conditions. Although the immune system normally removes damaged or abnormal cells from the body, some cancer cells are able to "hide" from the immune system.

Cancer is a genetic disease-that is, it is caused by changes to genes that control the way our cells function, especially how they grow and divide. The genetic changes that contribute to cancer tend to affect three main types of genes- proto-oncogenes, tumor suppressor genes and DNA repair genes. These changes are sometimes called "drivers" of cancer. Tumor markers are a kind of substances bio-synthesized and secreted by gene expression of tumor cells, or abnormally generated and fluctuating due to body responses toward tumors. They play an important role in early detection, diagnosis and prognosis of tumors.

Elabscience is passionate about cancer research; each of our ELISA Kits is created with careful consideration for the relevance and specificity of the targets as well as the necessity of validated applications. We are meticulous in the design and testing of our products before they reach you, and unwavering in our commitment to providing scientific technical support through the culmination of your project.



Invasion/Microenvironment

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0006	Human ANG (Angiogenin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0009	Human BAFF/CD257 (B-Cell Activating Factor) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0012	Human BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0013	Human BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0014	Human E-Cad (E-Cadherin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0017	Human CASP3 (Caspase 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0018	Human CTSB (Cathepsin B) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0020	Human MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0027	Human MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0048	Human IL-8 (Interleukin 8) ELISA Ki	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0050	Human IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0056	Human Cyt-C (Cytochrome C) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0060	Human EGFR (Epidermal Growth Factor Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0063	Human ES (Endostatin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0071	Human AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0079	Human G-CSF (Granulocyte Colony Stimulating Factor 3) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-H0084	Human HGF (Hepatocyte Growth Factor) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H0087	Human IGFBP-3 (Insulin Like Growth Factor Binding Protein 3) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0110	Human TGF- β 1 (Transforming Growth Factor Beta 1) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0111	Human VEGF-A (Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0116	Human LH (Luteinizing Hormone) ELISA Kit	0.94--60mIU/mL	0.56mIU/mL	Sandwich
E-EL-H0132	Human MMP-10 (Matrix Metalloproteinase 10) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0134	Human MMP-13 (Matrix Metalloproteinase 13) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0141	Human PRL (Prolactin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0151	Human IL-12 p40 (Interleukin 12 p40) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0175	Human HCG (Chorionic Gonadotropin) ELISA Kit	7.81--500mIU/mL	4.69mIU/mL	Sandwich
E-EL-H0179	Human FN (Fibronectin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0184	Human TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0237	Human AR (Amphiregulin) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0257	Human ACTN2 (Actinin Alpha 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0324	Human ANG (Angiostatin) ELISA Kit	31.25--2000ng/mL	18.75ng/mL	Sandwich
E-EL-H0422	Human ANXA5 (Annexin A5) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0425	Human IGF1R (Insulin Like Growth Factor 1 Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0435	Human IGF-2 (Insulin Like Growth Factor 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0483	Human bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	12.5--800pg/mL	7.5pg/mL	Sandwich
E-EL-H0513	Human MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0562	Human Bax (Bcl-2 Associated X Protein) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0570	Human bTC (Betacellulin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0669	Human CTSD (Cathepsin D) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-H0910	Human TP53 (Tumor Protein p53) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H0932	Human Smad1/MADH1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1051	Human GAL1 (Galectin 1) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1058	Human GAL7 (Galectin 7) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H1087	Human VEGFR1/FLT1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H1141	Human TSH (Thyroid Stimulating Hormone) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1298	Human SDC1 (Syndecan 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1347	Human OPN (Osteopontin) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1441	Human MMP-1 (Matrix Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1446	Human MMP-3 (Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1449	Human MMP-7 (Matrix Metalloproteinase 7) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1450	Human MMP-8 (Matrix Metalloproteinase 8) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1451	Human MMP-9 (Matrix Metalloproteinase 9) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1453	Human TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1454	Human TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1555	Human PGF (Placental Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1574	Human PTX 3/TSG-14 (Pentraxin 3, Long) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1584	Human Surv (Survivin) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1586	Human TGF- α (Transforming Growth Factor α) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1587	Human TGF- β 2 (Transforming Growth Factor β 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1590	Human TSP-2 (Thrombospondin-2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1593	Human TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1599	Human PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1600	Human VEGF-C (Vascular Endothelial Growth Factor C) ELISA Kit	109.38--7000pg/mL	65.63pg/mL	Sandwich
E-EL-H1606	Human EMMPRIN/CD147 (Extracellular Matrix Metalloproteinase Inducer) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1779	Human PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1855	Human KLK1 (Kallikrein 1) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1946	Human NE/ELA2 (Elastase 2, Neutrophil) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H2104	Human PAI1 (Plasminogen Activator Inhibitor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H2157	Human VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H2166	Human sVCAM-1/CD106 (soluble Vasuolar Cell Adhesion Molecule 1) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H2168	Human VWF (Von Willebrand Factor) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H2216	Human PF4 (Platelet Factor 4) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H2305	Human TACE (TNF α Converting Enzyme) ELISA Kit	62.5--4000pg/ml	37.5pg/ml	Sandwich
E-EL-H2312	Human TWEAK (Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H5200	Human LTF(Lactoferrin)ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H5439	Human INS (Insulin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H5587	Human VCAM-1/CD106 (Vasucular Cell Adhesion Molecule 1) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H5589	Human AITRL/TNFSF18(Activation-inducible TNF-related Ligand) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0004	Mouse BAFF/CD257 (B-cell Activating Factor) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0006	Mouse MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0013	Mouse MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0021	Mouse IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0023	Mouse Cyt C (Cytochrome C) ELISA Kit	0.39--25ng/mL	0.23ng/mL	Sandwich
E-EL-M0031	Mouse G-CSF (Granulocyte Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0035	Mouse IGFBP-3 (Insulin-like Growth Factor Binding Protein 3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0051	Mouse TGF- β 1 (Transforming Growth Factor β 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0057	Mouse LH (Luteinizing Hormone) ELISA Kit	0.47--30ng/mL	0.28ng/mL	Sandwich
E-EL-M0074	Mouse MMP-10 (Matrix Metalloproteinase 10) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0076	Mouse MMP-13 (Matrix Metalloproteinase 13) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0078	Mouse ANG (Angiogenin) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-M0083	Mouse PRL (Prolactin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0096	Mouse AGT (Angiotensinogen) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-M0097	Mouse ANG1 (Angiopietin 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0098	Mouse ANG2 (Angiopietin 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0099	Mouse ANXA5 (Annexin A5) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0170	Mouse bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0178	Mouse Bax (Bcl-2 Associated X Protein) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0181	Mouse BTC (Betacellulin) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0194	Mouse BMP-4(Bone Morphogenetic Protein 4)ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0196	Mouse BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0210	Mouse VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0211	Mouse E-Cad (E-Cadherin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0238	Mouse CASP3 (Caspase 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0249	Mouse CTSD (Cathepsin D) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0453	Mouse ES (Endostatin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0470	Mouse EPCAM (Epithelial Cell Adhesion Molecule) ELISA Kit	0.08--5ng/mL	0.05ng/mL	Sandwich
E-EL-M0506	Mouse FN (Fibronectin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M0527	Mouse GAL1 (Galectin 1) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-M0532	Mouse GAL7 (Galectin 7) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0626	Mouse MMP-3 (Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0627	Mouse MMP-9 (Matrix Metalloproteinase-9) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0629	Mouse OPN (Osteopontin) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-M0641	Mouse TIMP-1 (Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R1079	Rat VWF (Von Willebrand Factor) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0648	Mouse FLT1/VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0664	Mouse HGF (Hepatocyte Growth Factor) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0704	Mouse IGF-2 (Insulin Like Growth Factor 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0726	Mouse IL-12 (Interleukin 12) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0740	Mouse AR (Amphiregulin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0746	Mouse LTF/LF (Lactoferrin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0783	Mouse MMP-7 (Matrix Metalloproteinase 7) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0784	Mouse MMP-8 (Matrix Metalloproteinase 8) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0911	Mouse PGF (Placental Growth Factor) ELISA Kit	12.5--800pg/mL	7.5pg/mL	Sandwich
E-EL-M0914	Mouse PAI-1 (Plasminogen Activator Inhibitor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0926	Mouse PF4 (Platelet Factor 4) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M1062	Mouse Smad1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit	0.31--20ng/mL	0.13ng/mL	Sandwich
E-EL-M1064	Mouse ACTN2 (Actinin Alpha 2) ELISA Kit	0.31--20ng/mL	0.13ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M1084	Mouse TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1153	Mouse TSH (Thyroid Stimulating Hormone) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-M1164	Mouse TIMP-2 (Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M1165	Mouse TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M1190	Mouse TGF- α (Transforming Growth Factor α) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1226	Mouse PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M1230	Mouse VEGF-C (Vascular Endothelial cell Growth Factor C) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M1233	Mouse VCAM-1 (Vasculuar cell adhesion molecule 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M1247	Mouse vWF (von Willebrand Factor) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-M1274	Mouse PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M2413	Mouse AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M2423	Mouse CTSB (Cathepsin B) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2441	Mouse EGFR (Epidermal Growth Factor Receptor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R1057	Rat VEGF-C (Vascular Endothelial Cell Growth Factor C) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M2451	Mouse IL-12/p40 (Interleukin 12 p40) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M2460	Mouse SDC1 (Syndecan 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2614	Mouse INS(Insulin) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0006	Rat Cyt-C (Cytochrome C) ELISA Kit	0.39--25ng/mL	0.23ng/mL	Sandwich
E-EL-R0026	Rat LH (Luteinizing Hormone) ELISA Kit	1.56--100mIU/mL	0.94mIU/mL	Sandwich
E-EL-R0043	Rat MMP-10 (Matrix Metalloproteinase 10) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0045	Rat MMP-13 (Matrix Metalloproteinase 13) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0052	Rat PRL (Prolactin) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-R0063	Rat IL-12/p40 (Interleukin 12 p40) ELISA Kit	11.72--750pg/mL	7.03pg/mL	Sandwich
E-EL-R0064	Rat IL-12 (Interleukin 12) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0084	Rat TGF- β 1(Transforming Growth Factor β 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0091	Rat bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0098	Rat BAX (BCL-2 Associated X Protein) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0101	Rat BTC (Betacellulin) ELISA Kit	9.38--600pg/mL	5.63pg/mL	Sandwich
E-EL-R0116	Rat BMP-4 (Bone Morphogenetic Protein 4) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0119	Rat BMP-7 (Bone Morphogenetic Protein 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0130	Rat VE-Cadherin (Vascular Endothelial Cadherin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0160	Rat CASP3 (Caspase 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0171	Rat CTSD (Cathepsin D) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0347	Rat E-Cad (E-Cadherin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0364	Rat ES (Endostatin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0397	Rat GAL1 (Galectin 1) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-R0401	Rat GAL7 (Galectin 7) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0454	Rat G-CSF (Granulocyte Colony Stimulating Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0493	Rat HBEGF (Heparin Binding Epidermal Growth Factor Like Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0496	Rat HGF (Hepatocyte Growth Factor) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0530	Rat IGF-2 (Insulin Like Growth Factor 2) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-R0533	Rat IGFBP-3 (Insulin-Like Growth Factor Binding Protein 3) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0540	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0546	Rat IP-10/CXCL10 (Interferon Gamma Induced Protein 10kDa) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0574	Rat AR (Amphiregulin) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0577	Rat LTF/LF (Lactoferrin) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0578	Rat FN (Fibronectin) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-R0605	Rat MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0619	Rat MMP-3 (Matrix Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0622	Rat MMP-7 (Matrix Metalloproteinase 7) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0623	Rat MMP-8 (Matrix Metalloproteinase 8) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0624	Rat MMP-9(Matrix Metalloproteinase 9)ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-R0625	Rat ANG (Angiogenin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0633	Rat MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0637	Rat Smad1 (Mothers Against Decapentaplegic Homolog 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0702	Rat OPN (Osteopontin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-R0742	Rat PGF (Placental Growth Factor) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-R0744	Rat PAI-1 (Plasminogen Activator Inhibitor 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0759	Rat PF-4 (Platelet Factor-4) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0865	Rat TACE (TNF α Converting Enzyme) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0885	Rat ACTN2 (Actinin Alpha 2) ELISA Kit	1.0--25ng/mL	0.1ng/mL	Sandwich
E-EL-R0911	Rat FLT1/VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0976	Rat TSH (Thyroid Stimulating Hormone) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-R0985	Rat TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0986	Rat TIMP-3 (Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0996	Rat SDC1 (Syndecan 1) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-R1012	Rat TGF- α (Transforming Growth Factor Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1015	Rat TGF- β 2 (Transforming Growth Factor Beta 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1041	Rat TRAIL/TNFSF10 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1053	Rat PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R1055	Rat PLAUR/uPAR (Plasminogen Activator, Urokinase Receptor) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R1080	Rat XIAP (X-Linked Inhibitor Of Apoptosis Protein) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2386	Rat EGFR (Epidermal Growth Factor Receptor) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R2391	Rat ANXA5 (Annexin A5) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R2396	Rat BAFF/CD257 (B-Cell Activating Factor) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R2407	Rat AFGF/FGF1 (Acidic Fibroblast Growth Factor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R2414	Rat CTSB (Cathepsin B) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2466	Rat INS (Insulin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich

Tumor Immunology

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0032	Human CTACK (Cutaneous T-cell Attracting Chemokine) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0075	Human Flt3L (FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0088	Human IL-1 α (Interleukin 1 Alpha) ELISA Kit	3.91--250pg/mL	2.34pg/mL	Sandwich
E-EL-H0091	Human fPSA (Free Prostate Specific Antigen) ELISA Kit	0.31--20ng/m	0.19ng/mL	Sandwich
E-EL-H0099	Human IL-2 (Interleukin 2) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H0100	Human IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0101	Human IL-4 (Interleukin 4) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0102	Human IL-6 (Interleukin 6) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0103	Human IL-10 (Interleukin 10) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0104	Human IL-13(Interleukin 13)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0108	Human IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0109	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0149	Human IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0180	Human IL-9 (Interleukin 9) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0191	Human IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0222	Human IL-15 (Interleukin 15) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0235	Human IL-16 (Interleukin 16) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0643	Human CEA (Carcinoembryonic Antigen) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-H0648	Human IL-7 (Interleukin 7) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H1647	Human IL-12 p35 (Interleukin 12 p35) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H2338	Human IL-27 (Interleukin 27) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H2404	Human MIA1 (Melanoma Inhibitory Activity Protein 1) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H2450	Human IL-21 (Interleukin 21) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0030	Mouse Flt3L (FMS Like Tyrosine Kinase 3 Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0036	Mouse IL-1 α (Interleukin 1 Alpha)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0037	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-M0042	Mouse IL-2 (Interleukin 2) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0043	Mouse IL-4 (Interleukin 4) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0044	Mouse IL-6 (Interleukin 6) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0046	Mouse IL-10 (Interleukin 10) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0048	Mouse IFN- γ (Interferon Gamma) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0049	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0232	Mouse CEA (Carcinoembryonic Antign) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0371	Mouse CTACK (Cutaneous T-Cell Attracting Chemokine) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0721	Mouse IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0722	Mouse IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0723	Mouse IL-7 (Interleukin 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0724	Mouse IL-9 (Interleukin 9) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0727	Mouse IL-13 (Interleukin 13) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0728	Mouse IL-15 (Interleukin 15) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0729	Mouse IL-16 (Interleukin 16) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0732	Mouse IL-27 (Interleukin 27) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0961	Mouse PSA (Prostate Specific Antigen) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M1273	Mouse IL-21 (Interleukin 21) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M2450	Mouse IL-12 p35 (Interleukin 12 p35) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0009	Rat IFN- γ (Interferon Gamma) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0011	Rat IL-1 α (Interleukin 1 Alpha) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0012	Rat IL-1 β (Interleukin 1 Beta) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0013	Rat IL-2 (Interleukin 2) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0014	Rat IL-4 (Interleukin 4) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0015	Rat IL-6 (Interleukin 6) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0016	Rat IL-10 (Interleukin 10) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0019	Rat TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0150	Rat CEA (Carcinoembryonic Antigen) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0285	Rat CTACK (Cutaneous T-Cell Attracting Chemokine) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0556	Rat IL-3 (Interleukin 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0558	Rat IL-5 (Interleukin 5) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0559	Rat IL-7 (Interleukin 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0561	Rat IL-9 (Interleukin 9) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0563	Rat IL-13 (Interleukin 13) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0564	Rat IL-15 (Interleukin 15) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0565	Rat IL-16 (Interleukin 16) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0568	Rat IL-21 (Interleukin 21) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0769	Rat FXI (Coagulation Factor XI) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0796	Rat PSA (Prostate Specific Antigen) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2442	Rat IL-12 p35 (Interleukin 12 p35) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich

Growth Factors

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0074	Human FGF21 (Fibroblast Growth Factor 21) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0086	Human IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0092	Human FGF7 (Fibroblast Growth Factor 7) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0217	Human TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-H0442	Human IGFBP-1 (Insulin Like Growth Factor Binding Protein 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0443	Human IGFBP-5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0444	Human IGFBP-4 (Insulin Like Growth Factor Binding Protein 4) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0446	Human IGFBP-2 (Insulin Like Growth Factor Binding Protein 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0450	Human IGFBP-6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0552	Human FGF4 (Fibroblast Growth Factor 4) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0553	Human FGF6 (Fibroblast Growth Factor 6) ELISA Kit	93.75--6000pg/mL	56.25pg/mL	Sandwich
E-EL-H0748	Human CD30 (Cluster of Differentiation 30) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1117	Human FGF9 (Fibroblast Growth Factor 9) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1603	Human VEGFR-2/KDR (Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-H1629	Human EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-H1643	Human sTNF-1 (Soluble Tumor Necrosis Factor 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2214	Human PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0010	Mouse TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0029	Mouse FGF21(Fibroblast Growth Factor 21) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0034	Mouse IGF-1 (Insulin-like growth factor 1) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0171	Mouse FGF4 (Fibroblast Growth Factor 4) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0172	Mouse FGF6 (Fibroblast Growth Factor 6) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0285	Mouse CG (Chorionic Gonadotrophin) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-M0296	Mouse CD30 (Cluster of differentiation 30) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-M0499	Mouse FGF7 (Fibroblast Growth Factor 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0500	Mouse FGF9 (Fibroblast Growth Factor 9) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0649	Mouse VEGFR-2/KDR (Vasculoar Endothelial Growth Factor Receptor 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0705	Mouse IGFBP-1 (Insulin-Like Growth Factor Binding Protein 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0706	Mouse IGFBP-2 (Insulin-Like Growth Factor Binding Protein 2) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0707	Mouse PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-M0708	Mouse IGFBP-6 (Insulin-Like Growth Factor Binding Protein 6) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M1211	Mouse TRAIL R1 (Tumor Necrosis Factor-related Apoptosis-inducing Ligand 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2461	Mouse IGFBP5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-M2666	Mouse EGFR2 (Epidermal Growth Factor Receptor 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0010	Rat IGF-1 (Insulin-Like Growth Factor 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0088	Rat FGF4 (Fibroblast Growth Factor 4) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0089	Rat FGF6 (Fibroblast Growth Factor 6) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0090	Rat FGF9 (Fibroblast Growth Factor 9) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0216	Rat CD30 (Cluster of differentiation 30) ELISA Kit	0.16--100ng/mL	0.94ng/mL	Sandwich
E-EL-R0307	Rat FGF7 (Fibroblast Growth Factor 7) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R0421	Rat GTR (Glucocorticoid Induced Tumor Necrosis Factor Receptor) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-R0526	Rat IGFBP5 (Insulin Like Growth Factor Binding Protein 5) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0527	Rat IGFBP6 (Insulin Like Growth Factor Binding Protein 6) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0531	Rat IGFBP-1 (Insulin-Like Growth Factor Binding Protein 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0532	Rat IGFBP-2 (Insulin-Like Growth Factor Binding Protein 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0534	Rat IGFBP-4 (Insulin-Like Growth Factor Binding Protein 4) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-R1038	Rat TRAIL R1 (Tumor Necrosis Factor Related Apoptosis Inducing Ligand 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R1060	Rat VEGFR-2/KDR (Vascuoar Endothelial Growth Factor Receptor 2) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R1091	Rat TNFRSF1A (Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0084	Rat TGF- β 1(Transforming Growth Factor β 1) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0091	Rat bFGF/FGF2 (Basic Fibroblast Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R2408	Rat FGF21(Fibroblast Growth Factor 21) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich

Cancer Metabolism

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-0085	PGF (Prostaglandin F) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Competitive
E-EL-H0002	Human ACE (Angiotensin I Converting Enzyme) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0004	Human ADP/Acrp30 (Adiponectin) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H0007	Human ANG1 (Angiopoietin 1) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H0058	Human DPP4 (Dipeptidyl Peptidase IV) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0064	Human ET-1 (Endothelin 1) ELISA Kit	1.25--80pg/mL	0.75pg/mL	Sandwich
E-EL-H0066	Human EPO (Erythropoietin) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0093	Human LEPR (Leptin Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0142	Human MYO (Myoglobin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H0464	Human ApoB (Apolipoprotein B) ELISA Kit	78.13--5000ng/mL	46.88ng/mL	Sandwich
E-EL-H0467	Human ApoC3 (Apolipoprotein C3) ELISA Kit	7.81--500ng/mL	4.69ng/mL	Sandwich
E-EL-H0965	Human α 1-AT (Alpha 1-Antitrypsin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H1113	Human SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Competitive
E-EL-H1575	Human PDGF-AA (Platelet Derived Growth Factor-AA) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1817	Human GPI (Glucose 6 Phosphate Isomerase) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1818	Human GCK (Glucokinase) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-H1964	Human Anti-MPO (Anti-Myeloperoxidase Antibody) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-H2172	Human HO1 (Heme Oxygenase 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H2420	Human PARP4 (Poly ADP Ribose Polymerase 4) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H5408	Human CAT(Catalase) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0002	Mouse ADP/Acrp30 (Adiponectin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0027	Mouse EPO (Erythropoietin) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0084	Mouse MYO (Myoglobin) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-M0132	Mouse ApoB (Apolipoprotein B) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0457	Mouse ET-1 (Endothelin 1) ELISA Kit	7.81--500ng/mL	4.69ng/mL	Sandwich
E-EL-M0559	Mouse GCK (Glucokinase) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0565	Mouse GPI (Glucose 6 Phosphate Isomerase) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0630	Mouse PDGF-AA (Platelet-Derived Growth Factor-AA) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M0651	Mouse HO-1 (Heme Oxygenase 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0749	Mouse LEPR (Leptin Receptor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M2398	Mouse SOD1 (Superoxide Dismutase 1, Soluble) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Competitive
E-EL-M2428	Mouse α 1-AT (Alpha 1-Antitrypsin) ELISA Kit	3.13--200ng/mL	1.88ng/mL	Sandwich
E-EL-M2440	Mouse DPP4 (Dipeptidyl Peptidase IV) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M2463	Mouse ApoC3 (Apolipoprotein C3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M2605	Mouse CAT (Catalase) ELISA kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M2632	Mouse MPO(Myeloperoxidase)ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0007	Rat EPO (Erythropoietin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0053	Rat MYO (Myoglobin) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-R0167	Rat ET-1 (Endothelin 1) ELISA Kit	7.81--500ng/mL	4.69ng/mL	Sandwich
E-EL-R0329	Rat ADP/Acrp30 (Adiponectin) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-R0337	Rat DPP IV (Dipeptidyl Peptidase IV) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-R0426	Rat GCK (Glucokinase) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-R0488	Rat HO-1 (Heme Oxygenase 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0626	Rat ANG1 (Angiopietin 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R1095	Rat LEPR (Leptin Receptor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R1218	Rat ApoB (Apolipoprotein B) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R1424	Rat SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Competitive
E-EL-R2379	Rat ApoC3 (Apolipoprotein C3) ELISA Kit	6.25--400ng/ml	3.75ng/ml	Sandwich
E-EL-R2401	Rat ACE (Angiotensin I Converting Enzyme) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-R2420	Rat α 1-AT (Alpha 1-Antitrypsin) ELISA Kit	6.25--400ng/ml	3.75ng/ml	Sandwich
E-EL-R2456	Rat CAT(Catalase) ELISA kit	15.63--1000pg/mL	9.38pg/mL	Sandwich

Signal Transduction

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0082	Human sgp130 (soluble Glucoprotein 130) ELISA Kit	0.25--16ng/mL	0.15ng/mL	Sandwich
E-EL-H0493	Human AR (Androgen Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H0930	Human SHH (Hedgehog Homolog, Sonic) ELISA Kit	78.13--5000ng/mL	46.88ng/mL	Sandwich
E-EL-H1133	Human ER α (Estrogen Receptor Alpha) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1174	Human PKB (Protein Kinase B) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1324	Human Smad3/MADH3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-H2130	Human STAT1 (Signal Transducer and Activator of Transcription 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H2132	Human STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H2271	Human PGR (Progesterone Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0476	Mouse ER α (Estrogen Receptor Alpha) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-M0509	Mouse FAK (Focal Adhesion Kinase) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-M0560	Mouse gp130 (Glucoprotein 130) ELISA Kit	125--8000pg/mL	75pg/mL	Sandwich
E-EL-M0623	Mouse SHH (Hedgehog Homolog, Sonic) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0944	Mouse PGR (Progesterone receptor) ELISA Kit	1.0--25ng/mL	0.1ng/mL	Sandwich
E-EL-M0969	Mouse PIAS1 (Protein Inhibitor Of Activated STAT 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0972	Mouse PKB (Protein Kinase B) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-M1057	Mouse STAT3 (Signal Transducer And Activator Of Transcription 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2648	Mouse sgp130(Soluble Glucoprotein 130)ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0211	Rat JNK (C-Jun N-Terminal Kinases) ELISA Kit	7.81--500ng/mL	4.69ng/mL	Sandwich
E-EL-R0484	Rat SHH (Hedgehog Homolog, Sonic) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Competitive
E-EL-R0713	Rat Smad3 (Mothers Against Decapentaplegic Homolog 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0776	Rat PGR (Progesterone Receptor) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0807	Rat PKB (Protein Kinase B) ELISA Kit	3.13--200U/ml	1.88U/ml	Sandwich
E-EL-R1226	Rat STAT3 (Signal Transducer and Activator of Transcription 3) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich



Oncoproteins/Suppressors

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0068	Human FASL/TNFSF6 (Factor Related Apoptosis Ligand) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0096	Human NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H1161	Human PRDX1 (Peroxisredoxin 1) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H2367	Human DKK3 (Dickkopf Related Protein 3) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H2473	Human CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0028	Mouse FASL (Factor Related Apoptosis Ligand) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0118	Mouse CDKN1A (Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit	7.81-500pg/mL	1.95pg/mL	Sandwich
E-EL-M0443	Mouse DKK3 (Dickkopf Related Protein 3) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0828	Mouse NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-M1215	Mouse TP53 (Tumor Protein 53) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R0310	Rat CDKN1A (Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R0645	Rat PTEN (Phosphatase and Tensin Homolog) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R0662	Rat NGAL (Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-R1042	Rat TP53 (Tumor Protein 53) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-R2373	Rat FASL (Factor Related Apoptosis Ligand) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R2452	Rat DKK3 (Dickkopf Related Protein 3) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich

Tumor Biomarkers

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0062	Human ENG (Endoglin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0070	Human αFP (Alpha-Fetoprotein) ELISA Kit	1.56--100ng/mL	0.94ng/mL	Sandwich
E-EL-H0112	Human EREG (Epregulon) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H0636	Human CA125 (Carbohydrate Antigen 125) ELISA Kit	3.13--200U/mL	1.88U/mL	Sandwich
E-EL-H0637	Human CA19-9 (Carbohydrate antigen19-9) ELISA Kit	6.25--400U/mL	3.75U/mL	Sandwich
E-EL-H1094	Human VIM (Vimentin) ELISA Kit	7.81-500pg/mL	4.69ng/mL	Sandwich
E-EL-H1260	Human NNE (Non-Neuronal Enolase) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H1270	Human PARC/CCL18 (Pulmonary Activation Regulated Chemokine) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H5309	Human sENG (Soluble Endoglin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M0870	Mouse CA125 (Carbohydrate Antigen 125) ELISA Kit	3.13--200U/mL	1.88U/mL	Sandwich
E-EL-M0991	Mouse PARC (Pulmonary Activation Regulated Chemokine) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M1072	Mouse ENG (Endoglin) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-M2405	Mouse αFP (Alpha-Fetoprotein) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2427	Mouse EREG (Epregulin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0153	Rat αFP(Alpha-Fetoprotein)ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0707	Rat CA125 (Carbohydrate Antigen 125) ELISA Kit	7.81--500U/mL	4.69U/mL	Sandwich
E-EL-R2403	Rat NNE (Non-Neuronal Enolase) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R2419	Rat EREG (Epregulin) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R2434	Rat ENG (Endoglin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-R2539	Rat CA19-9 (Carbohydrate antigen19-9) ELISA Kit	1.56--100U/mL	0.94U/mL	Sandwich

Cell Death

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0034	Human FcεRII/CD23 (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0059	Human EGF (Epidermal Growth Factor) ELISA Kit	3.91--250pg/mL	2.34pg/mL	Sandwich
E-EL-H0094	Human LIF (Leukemia Inhibitory Factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H0479	Human AIF (Apoptosis Inducing Factor) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H0671	Human CTSL (Cathepsin L) ELISA Kit	62.5--4000pg/mL	37.5pg/mL	Sandwich
E-EL-H1205	Human NGF (Nerve Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-H1576	Human PDGF-AB (Platelet Derived Growth Factor AB) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-H1634	Human FXN (Frataxin) ELISA Kit	23.44--1500pg/mL	14.06pg/mL	Sandwich
E-EL-H1655	Human mTOR (Mammalian Target of Rapamycin) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-H2086	Human ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-H2087	Human ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit	46.88--3000pg/mL	28.13pg/mL	Sandwich
E-EL-H5022	Human IL-11 (Interleukin 11) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0025	Mouse EGF (Epidermal growth factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0040	Mouse LIF (Leukemia Inhibitory Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0141	Mouse AIF (Apoptosis Inducing Factor) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich



Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-M0251	Mouse CTSL (Cathepsin L) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0480	Mouse ERK1 (Extracellular Signal Regulated Kinase 1) ELISA Kit	78.13--5000pg/mL	46.88pg/mL	Sandwich
E-EL-M0481	Mouse ERK2 (Extracellular Signal Regulated Kinase 2) ELISA Kit	5.0--100ng/mL	1.0ng/mL	Sandwich
E-EL-M0631	Mouse PDGF-AB (Platelet-Derived Growth Factor AB) ELISA Kit	78.12-5000pg/mL	32pg/mL	Sandwich
E-EL-M0725	Mouse IL-11 (Interleukin 11) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-M0815	Mouse NGF (Nerve growth factor) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-M0876	Mouse MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-M2432	Mouse FcεRII/CD23 (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R0173	Rat CTSL (Cathepsin L) ELISA Kit	0.63--40ng/mL	0.38ng/mL	Sandwich
E-EL-R0369	Rat EGF (Epidermal Growth Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0372	Rat ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit	0.31--20ng/mL	0.19ng/mL	Sandwich
E-EL-R0536	Rat PDGF-AB (Platelet-Derived Growth Factor AB) ELISA Kit	7.81--500pg/mL	4.69pg/mL	Sandwich
E-EL-R0562	Rat IL-11 (Interleukin 11) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0585	Rat LIF (Leukemia Inhibitory Factor) ELISA Kit	15.63--1000pg/mL	9.38pg/mL	Sandwich
E-EL-R0652	Rat NGF (Nerve Growth Factor) ELISA Kit	39.06--2500pg/mL	23.44pg/mL	Sandwich
E-EL-R0711	Rat MAPK14 (Mitogen Activated Protein Kinase 14) ELISA Kit	1.25--80ng/mL	0.75ng/mL	Sandwich
E-EL-R2422	Rat FcεRII/CD23 (Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	31.25--2000pg/mL	18.75pg/mL	Sandwich
E-EL-R2528	Rat ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich
E-EL-R2529	Rat ERK1/2(Extracellular Signal Regulated Kinase 1/2) ELISA Kit	0.78--50ng/mL	0.47ng/mL	Sandwich

Cancer Biomarker Panels

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H0616	Human MUC1 (Mucin 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich
E-EL-M2604	Mouse MUC1 (Mucin 1) ELISA Kit	0.16--10ng/mL	0.09ng/mL	Sandwich



Featured Publications

Cat. No.	Targets	Species	Publications
E-EL-H2050	pepsin	Human	Kang H S, Kim J K, Lee H J, et al. LRP1-dependent pepsin clearance induced by 2'-hydroxycinnamaldehyde attenuates breast cancer cell invasion[J]. The international journal of biochemistry & cell biology, 2014, 53: 15-23.
E-EL-M0253	caveolin-1	Mouse	Xu W, Nan S, Bai W, et al. Caveolin-1 as a biomarker to predict therapeutic effect of low-frequency ultrasound combined with SonoVue on prostate cancer in nude mice model[J]. Cancer Biomarkers, 2014, 14(4): 279-286.
E-EL-H0375	Soluble fatty acid synthase	Human	Qin W, Raatz S, Zhang K K, et al. Polyunsaturated Fatty Acid Content May Be Increased in the Milk of Women with Pregnancy-Associated Breast Cancer[J]. Journal of Human Lactation, 2014: 0890334414539576.
E-EL-H0190	C5a ELISA kit	Human	Wu Y, Wang Y, Qin F, et al. CD55 limits sensitivity to complement-dependent cytolysis triggered by heterologous expression of α -gal xenoantigen in colon tumor cells[J]. American Journal of Physiology-Gastrointestinal and Liver Physiology, 2014, 306(12): G1056-G1064.
E-EL-M0372	β -CTX	Mouse	Yan Z, Jin S, Wei Z, et al. Discoidin domain receptor 2 facilitates prostate cancer bone metastasis via regulating parathyroid hormone-related protein[J]. Biochimica et Biophysica Acta (BBA)-Molecular Basis of Disease, 2014, 1842(9): 1350-1363.
E-EL-H0091 E-EL-H0102 E-EL-H5413	PSA IL-6 PSMA	Human	Liu J, Lu C Y, Zhou H, et al. Flexible gold electrode array for multiplexed immunoelectrochemical measurement of three protein biomarkers for prostate cancer[J]. ACS applied materials & interfaces, 2014, 6(22): 20137-20143.
E-EL-H1274	Human SEMA3A	Human	Wang Z, Chen J, Zhang W, et al. Axon guidance molecule semaphorin3A is a novel tumor suppressor in head and neck squamous cell carcinoma[J]. Oncotarget, 2016, 7(5): 6048.
E-EL-M0890 E-EL-M2664	Mouse PRF1 (Perforin 1) Mouse LAMP1	Mouse	Zhang H, Li N, Zhang J, et al. The Influence of miR-34a expression on stemness and cytotoxic susceptibility of breast cancer stem cells[J]. Cancer biology & therapy, 2016: 1-11.

Please click the link www.elabscience.com for more publications



International Quality and Management
Certification ISO9001:2008



Elabscience®

